

Spring 2010 Microscope Deals

all display units **40%** off list



AE31

List: \$3,495
Sale: **2,095**



B5

List:	Binocular \$1,595	Trinocular \$1,895
Sale:	955	1,135

Research Grade Inverted Microscope

AE31 Trinocular (80/20) Inverted Microscope, WF10x/22 Eyepieces, Plan Achromat PL4x Plan Achromat Phase 10x and 20x objectives, ELWD condenser, Phase slider, PH1, PH3 and phase centering telescope, plain stage, 6V/30W Quartz illumination

Professional Series Compound Microscope

Infinity Optical System Parfocal Distance 45mm, Siedentopf binocular, quituple nosepiece, Abbe condenser N.A. 1.25, and 30.0 mm widefield eyeieces,



BA450

List: \$2,775
Sale: **1,665**



BA200

Teaching scopes



BA300

List: \$995
Sale: **597**

List: \$1,450
Sale: **870**

Research & Clinical Biological Microscope

Color Corrected Infinity Optical System, CCIS optics allow new accessories and functions to be incorporated into the BA400/450 clinical microscope. Motic infinity objectives have succeeded in achieving longer working distances with higher numerical aperture while significantly improving optical performance and versatility.



SMZ-140

Zoom Stereomicroscope

List: \$850
Sale: **510**



K-700L

Stereomicroscope

List: \$1,400
Sale: **840**



SCIMETRICS INC.

LABORATORY INSTRUMENTATION "Equipping Labs for Life"

19407 Park Row #102 Houston, TX 77084
www.scimetricsinc.com scimetrics@scimetricsinc.com
(800) 231-2065 FAX: (281) 565-0570