



Hermle Z400 Universal Centrifuge

Maximum speed 13,500rpm/17,310 x g. Automatic soft stop braking. Microprocessor control with g-force conversion program. Multiple swing out and angle rotor options. Features maintenance-free induction drives that provide exceptionally quiet and smooth operation.

This is a display unit - never used.

List Price: \$4,300

Sale Price: \$3,200

Please specify rotors and inserts needed. All are available at 20% off list.



C2400

C1301-B

Maximum speed 13,300rpm/16,000 x g.
Complete with easy access, 24 x 1.5/2.0 ml rotor
Set and view speed in rpm or g-force
Digital set and display of parameters
Quiet and cool running
Display unit available at **\$1,395** w/o iPod

Spectrafuge Mini Centrifuge



List Price: \$305

Sale Price: \$218

Maximum speed/RCF: 6,000/2,000
Supplied with both standard microtube and strip tube rotors
Ideal for quick spin downs and microfiltration
Starts and stops in seconds
Compact design
One display unit with blue lid at: **\$195**

Call or e-mail Scimetrics, Inc. to order